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- Análisis FMEA (Failure Mode Effects Analysis).

RISK IDENTIFICATION		RISK ANALYSIS		RISK EVALUATION		TOTAL	RISK CONTROL STRATEGY	
PRODUCT & PROCESS PERSPECTIVE								
PRODUCT CHARACTERIZATION	POTENTIAL RISK	DETECTABILITY	SEVERITY	PROBABILITY	RISK PRIORITY NUMBER	RISK REDUCTION METHOD	FINAL RISK	COMMENTS
Product Toxicity OEL OEL <0,01µg/m3	Exposure during sampling/weighting	3	3	3	27	Use Isolator	Accepted	Exposure <0,01µg/m3
Number Batch per day - Value 2	Total Exposure Risk Time, per day.	2	2	2	8	Use Isolator	Accepted	
Batch Size - Value 0.5 to 7 kg	Bigger batch, bigger risk.	2	2	2	8	Use Isolator	Accepted	
Product Persistence	Product remaining floating in the air	3	3	3	27	Use Isolator	Accepted	
PRIMARY BARRIER: MATERIAL TRANSFER PERSPECTIVE								
RAW MATERIAL PACKAGING	POTENTIAL RISK	DETECTABILITY	SEVERITY	PROBABILITY	RISK PRIORITY NUMBER	RISK REDUCTION METHOD	FINAL RISK	COMMENTS
Drums	Product loses during opening and close the container. Different drum sizes.	2	3	2	12	Drum Docking System compatible with different drums sizes (see Picture XX)	Accepted	Raw material supplier have to provide double seal and container bag inside.
Cans with threaded cap.	Product loses during opening and close the container	2	3	2	12	Use Air Lock chamber big enough, and open the cans inside the Isolator.	Accepted	Only valid for small containers.
Bottles, cans to transfer the product to production.	Product loses during opening and close the container	3	3	2	18	Spill Entry Valve or Single use full transfer System with special connections.	Accepted	Useful for transfer the product to production.
PRIMARY BARRIER: HUMAN RESOURCES PERSPECTIVE								
REGULAR PRODUCTION	POTENTIAL RISK	DETECTABILITY	SEVERITY	PROBABILITY	RISK PRIORITY NUMBER	RISK REDUCTION METHOD	FINAL RISK	COMMENTS
Normal Operation	Mistakes due to lack of training.	3	3	3	27	Training in HAPI sampling and weighing, GMP, Containment. Use of Mix Ready Technology.	Accepted	
Docking Systems and RTP	Mistakes due to lack of training in transfer systems.	3	3	2	18	Training to use transfer systems, docking systems and RTP connectors.	Accepted	
Primary barrier use.	Mistakes due to lack of training in primary barriers and Containment requirements.	3	3	2	18	Training to use Isolators.	Accepted	
PPE use.	PPE selection, clothes, gloves, mask, cap. Risk during gowning and degowning during the entrance and exit from containment area.	2	3	2	12	Training about PPE in containment environment.	Accepted	
Exposure during the Degowning process.	After an incident where an operator has been exposed to an abnormally large amount of product. Risk during degowning.	3	3	2	18	Mist shower before degowning. Emergency management training.	Accepted	
PRIMARY BARRIER: CLEANING AND MAINTENANCE PERSPECTIVE								
MAINTENANCE	POTENTIAL RISK	DETECTABILITY	SEVERITY	PROBABILITY	RISK PRIORITY NUMBER	RISK REDUCTION METHOD	FINAL RISK	COMMENTS
Open Isolator	Exposure of maintenance staff	3	3	3	27	Emergency management training, and containment training.	Accepted	
Filter Change	Exposure to the powders in the filters.	3	3	2	18	Use BMD filters.	Accepted	
Isolator Ventilation System Monitoring	Loss tightness and containment, for any technical reason.	2	3	2	12	Monitoring: pressure inside isolator, inflatable gasket status. Periodically HEPA filters integrity	Accepted	Including an Alarm management system.
CLEANING PROCESS	POTENTIAL RISK	DETECTABILITY	SEVERITY	PROBABILITY	RISK PRIORITY NUMBER	RISK REDUCTION METHOD	FINAL RISK	COMMENTS
Regular Cleaning	Exposure of staff. Cross contamination.	3	3	3	27	Primary Barrier Cleanability, rounded corners, drainability. Cleaning Validation.	Accepted	
WIP & Spray Balls	Insufficient cleaning possibilities, manual process.	3	3	3	27	Wet cleaning, manual and semiautomatic cleaning process with spray balls. Cleaning Validation.	Accepted	
Waste Material	Exposure during waste removal.	3	3	2	18	Safe waste material transfer. Containes liner or disposable betabag with RTP.	Accepted	
Waste Water	Product loses through drain point, after WP.	3	3	3	27	Connect drain point to a controlled drain network.	Accepted	

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RISK PRIORITY NUMBER (RPN)	PRIORITIES FOR ACTION AND DEDICATION OF RESOURCES
19 < RPN < 27 : HIGH RISK	Priority and immediate action Maximum dedication of resources and information
9 < RPN < 18: MODERATE RISK	Secondary action, through an action plan Less dedication of resources and information
1 < RPN < 8. LOW RISK	Deferred action or not required Minimum dedication of resources and information

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1. ICH guideline Q9 on quality risk management.
2. ISPE Baseline Guide Volume 1- Active Pharmaceutical Ingredients.
3. ISPE Good Practice Guide - Assessing the particulate Containment Performance of Pharmaceutical Equipment (SMEPAC).
4. ISPE Containment Manual- D/A/CH Affiliate 2016.
5. Ader A.W., Farris J.P., Ku R.H., Occupational health, categorization and compound handling practice systems - roots, application and future. *Chemical Health and Safety*, 12(4), 20-26 (2005).
6. Farris J.P., Ader A.W., Ku R.H., History, implementation and evolution of the pharmaceutical hazard categorization and control system, *Chimica oggi - Chemistry Today*, 24(2), 5-10 (2006).
7. SafeBridge Consultants, Inc., Occupational Health Toxicity/Potency Categorization and Handling Practices - SafeBridge, Fifth Revision - January 2002, pp. 12.
8. Naumann B.D., Sargent E.V., Starkman B.S., et al., Performancebased exposure control limits for pharmaceutical active ingredients, *American Industrial Hygiene Association Journal*, 57(1), 33-42 (1996).
9. Naumann B.D. Control banding in the pharmaceutical industry: www.aioh.org.au/downloads/documents/ControlBandingBNaumann.pdf (last checked on Feb. 3rd 2014).
10. Manufacturing Highly Potent Drugs, reducing the risks, *Pharmaceutical Visions*, SafeBridge.
11. Markets and Markets, "High Potency APIs /HPAPI Market worth 26.84 Billion USD by 2023," Press Release, 04/18.

..... 20 - 21 SEPTIEMBRE 2023

PABELLÓN 14
IFEMA
MADRID

Cosméticaforu!!!

farmaforu!!!

Labforu!!!

..... CON LAS ÁREAS ESPECIALIZADAS

Health+Researchforu!!!

Biotechforu!!!

Innovaforu!!!

APIsforu!!!

Cannabisforu!!!

Nutraforu!!!

INFORMACIÓN

☎ +34 91 630 85 91 / +34 672 050 625 ✉ info@farmaforum.es 🌐 farmaforum.es